

## How to integrate transdermal-vitals-scanner library to android app

1. Add tvs.shared-1.0.0.jar library to project and add as local dependency (for example, in gradle something like `implementation(files("libs/com.tvs.shared-1.0.0.jar"))`)

2. Add these configuration to access camera to your AndroidManifest file in resource folder:

```
<uses-permission android:name="android.permission.CAMERA" />
<uses-feature android:name="android.hardware.camera" />
<uses-feature android:name="android.hardware.camera.flash" />
<uses-feature android:name="android.hardware.camera.autofocus" />
<uses-permission android:name="android.permission.FLASHLIGHT" />
<uses-permission android:name="android.permission.WAKE_LOCK" />
```

After that you need to check and request camera permissions if needed

```
private boolean checkPermission() {
    return ContextCompat.checkSelfPermission(this, Manifest.permission.CAMERA) == PackageManager.PERMISSION_GRANTED;
}

private void requestPermission() {
    ActivityCompat.requestPermissions(this, new String[]{Manifest.permission.CAMERA}, PERMISSION_REQUEST_CODE);
}
```

3. There is one method `processImage` in `vitals.*Processor` files which you should use in your android app. It accepts the byte array with image data, camera size (weight + height) and some anthropometric data. It returns `ProcessStatus` enum, with 5 statuses

```
RED_INTENSITY_NOT_ENOUGH("Not good red intensity to process. Should start again"),
MEASUREMENT_FAILED("Measurement Failed. Should Start again"),
IN_PROGRESS("Processing in progress"),
PROCESS_FINISHED("Processing finished"),
NEED_MORE_IMAGES("Need more images to process")
```

4. In your codebase with android components you should implement class extends `AppCompatActivity` which opens camera and process images in `onPreviewFrame` method. Example could be found below:

```
class VitalSignsProcess : AppCompatActivity() {
    private var preview: SurfaceView? = null
    private lateinit var vsp: VitalSignsProcessor

    private var previewHolder: SurfaceHolder? = null
    private var camera: Camera? = null
    private var wakeLock: PowerManager.WakeLock? = null
    private var mainToast: Toast? = null
    private var progBar: ProgressBar? = null
    var progP = 0
```

```

var inc = 0

override fun onCreate(savedInstanceState: Bundle?) {
    super.onCreate(savedInstanceState)
    setContentView(R.layout.activity_vital_signs_process)
    vsp = VitalSignsProcessor(
        User(
            height = 188.0,
            weight = 76.0,
            age = 25,
            gen = 1
        )
    )
    preview = findViewById(R.id.preview)
    previewHolder = preview!!.holder
    previewHolder!!.addCallback(surfaceCallback)
    previewHolder!!.setType(SurfaceHolder.SURFACE_TYPE_PUSH_BUFFERS)
    progBar = findViewById(R.id.VSPB)
    progBar!!.progress = 0
    val pm: PowerManager = getSystemService(Context.POWER_SERVICE)
as PowerManager
    wakeLock = pm.newWakeLock(PowerManager.FULL_WAKE_LOCK, "tvs:
DoNotDimScreen")
    }

    //getting frames data from the camera and start the measuring
process
    private val previewCallback: Camera.PreviewCallback = object :
Camera.PreviewCallback {

        override fun onPreviewFrame(data: ByteArray, cam: Camera) {
            val size = cam.parameters.previewSize

            val width = size.width
            val height = size.height

            when (vsp.processImage(data, width, height)) {
                ProcessStatus.RED_INTENSITY_NOT_ENOUGH -> {
                    inc = 0
                    progP = inc
                    progBar!!.progress = progP
                }
                ProcessStatus.MEASUREMENT_FAILED -> {
                    inc = 0
                    progP = inc
                    progBar!!.progress = progP
                    mainToast =
                        Toast.makeText(applicationContext, "Measurement
Failed", Toast.LENGTH_SHORT)
                    mainToast!!.show()
                }
            }
        }
    }
}

```

```

        }
        ProcessStatus.IN_PROGRESS -> {
            progP = inc++ / VITALS_PROCESS_DURATION
            progBar!!.progress = progP
        }
        ProcessStatus.PROCESS_FINISHED -> {
            //do what you need. Data can be accessed from vsp.*
fields (vsp.o2, vsp.Breath, vsp.Beats, vsp.SP, vsp.DP)
            finish()
        }
        ProcessStatus.NEED_MORE_IMAGES -> {
            //do nothing. Need to process more images
        }
    }
}

private val surfaceCallback: SurfaceHolder.Callback = object :
SurfaceHolder.Callback {
    override fun surfaceCreated(holder: SurfaceHolder) {
        try {
            camera!!.setPreviewDisplay(holder)
            camera!!.setPreviewCallback(callback)
        } catch (t: Throwable) {
        }
    }

    override fun surfaceChanged(holder: SurfaceHolder, format: Int,
width: Int, height: Int) {
        val parameters = camera!!.parameters
        parameters.flashMode = Camera.Parameters.FLASH_MODE_TORCH
        val size = getSmallestPreviewSize(width, height, parameters)
        if (size != null) {
            parameters.setPreviewSize(size.width, size.height)
        }
        camera!!.parameters = parameters
        camera!!.startPreview()
    }

    override fun surfaceDestroyed(holder: SurfaceHolder) {}
}

private fun getSmallestPreviewSize(
    width: Int,
    height: Int,
    parameters: Camera.Parameters
): Camera.Size? {
    var result: Camera.Size? = null
    for (size in parameters.supportedPreviewSizes) {
        if (size.width <= width && size.height <= height) {

```

```
        if (result == null) {
            result = size
        } else {
            val resultArea = result.width * result.height
            val newArea = size.width * size.height
            if (newArea < resultArea) result = size
        }
    }
}
return result
}
}
```

5. Results could be accessed from the \*Processor object fields (for example, vsp.o2, vsp.Breath, vsp.Beats, vsp.SP, vsp.DP)